

Automatically sends results via Bluetooth to our LifeGuard SmartHome and SmartKiosk systems



Features

- Measurements can be sent automatically via Bluetooth to our LifeGuard SmartHome and LifeGuard SmartKiosk TeleHealth systems for monitoring by a Nurse in our Emergency Response Centre
- Irregular heart beat detection
- Automated averaging function
- Comfort inflation technology
- Data transmission
- 400 memory sets
- 4 users selection

Specifications

- Accuracy of pressure: ± 3 mmHg or $\pm 2\%$ of reading
- Accuracy of HHeart Rate: $\pm 4\%$ of reading
- Measurement Unit: Systolic: 50-250 mmHg
Diastolic: 30-180 mmHg
- Measurement Range: Pulse Rate: 40-199 / minute
- Memory Capacity: 400 sets
- Data output: Bluetooth
- Power source: 4 x AA + 6V DC adaptor
- Dimension: 141 (L) x 121 (W) x 72 (H) mm
- Weight: 350g (without battery)
- Operating condition: 10°C - 40°C, below 85% R.H.

Specifications subject to change at discretion of INS LifeGuard or the manufacturer, TaiDoc Technology Corp.

Rev. 03/06/14



www.theinsgroup.com.au

1 800 636 040